


# COSMETIC GYNEGOLOGY quipment



# 1470nm / 15 W Diode Laser with LVR and SUI Probe For Cosmetic Gynecology



**SUI** (Stress Urinary Incontinence) Handpiece

#### THE BENEFITS OF LASER TREATMENT

#### A Simple Non Surgical Process that result in

- Increased Collagen
- Improved Muscle Tone
- Improved Sensation
- Reduced Vaginal Dryness
- Enhanced Orgasms
- OPD process without anesthesia
- No side effects
- Effective and painless

#### **APPLICATIONS**

- Vaginal rejuvenation
- GSM-Genitourinary syndrome of menopause
- Vaginal atrophy
- Vaginal laxity
- Diseases associated to postpartum alteration
- HPV
- Dryness
- Itching
- Vaginismus

#### SPECIFICATION:

- Laser type
- Wavelength
- MAX Power
- Aiming Beam

Pulse Time

- Operating Mode
- Beam delivery
- Optical Fiber Connector
- Emmition Initiation.
- Controller
- Display

- Diode, Semiconductor
- : 1470nm
- : 15 Watts
- 635nm/Pulse, Single, Group
- Continuous Pulse, Single, Group
- 0.05ms-1000ms
- : Fiber Optic with SMA 905 Connector
- : SMA 905, accepts optical fibers having
  - a core diameter from 200µm up to 600pm,
- : Footswitch
- : Microprocessor
- : 10.1" IPS with touch panel
- Internal, air and thermoelectric cooling Cooling System

- Power Supply of the Laser
- Power Supply of AC Adapter
- AC Adapter
- Laser Dimensions
- Laser weight
- Laser Case Dimensions
- Weight of Laser with Cases
- Class of Medical Device
- Laser Safety Class
- Electric Safety Class
- Environmental Conditions During Work

- : DC 24V/8.33A from the Separate AC
- : Single phase 100-240 VAC: 50-60 Hz;
- : DC 24v/8,33A Medical device class EN60601-1 Class |
- : 27cm\* 11cm\*19 cm length\*width\*height
- : 2.75kg
- : 56cm × 35cm × 23cm
- : 09 kg
- : IIB
- : 4
- : I type B
- : From +10 to 240C degree, relative humidity from 30% up to 60%



#### BENEFIT:

- Minimally Invasive, No Incision (no cutting, no bleeding, no ablation)
- Non-anesthesia, Painless
- Don't need any special care before & after procedure
- No Consumables
- Easy and Fast procedure
   High rate of success with only one treatment

#### SPECIFICATIONS:

Laser Type : CO2 laser


Wave Length : 10.6 um

Watt : 1 - 50W

Electric Requirements : 220 – 240 VAC, 50/60Hz

Dimension : 356(W) X 428(D) X 1260(H) mm

Weight : 45 Kg





# DIFFERENT MODES:

Surgical Mode : Surgery, Nevus, Wart, Tumor removal

Zoom Mode : Similar to surgical Application, with bigger spot size

Fractional Mode : Aesthetic, Acne scar, Pore, Wrinkle, Laser resurfacing, Pudenda whitening

Skin/Pudenda texture improvement

Beauty Pulse Mode : Vaginal rejuvenation, Vaginal elasticity Improvement, Colpoxerosis and SUI

improvement



# RF (Radio Frequency)

- Vaginal Rejuvination
- Majora/ Minora Devlopment
- Skin Lustering



# SPECIFICATION:

Input Power :  $100-120/200-240V \pm 10\%$ , 50/60~Hz Output Power : 20~watts (maximum)  $50-300~\Omega~0-10~volts$ ,

in Stimulate Mod

Maximum Output Voltage : 80V, RMS

Operating Frequency: 460 kHz (± 5kHz)

Set Temperature Range : 35-47°C; this range varies for each Power

Delivery Mode

Dimensions : 5.75" H X 12.5" W X 13.0" D

(146mm X 318mm X 330mm)

Weight : 12.6 lbs (5.7kg)

# HIFU (HIGH INTENSITY FOCUSED ULTRASOUND)



- 1) Sonic Vaginal rejuvenation
- 2) Sonic Vaginal tightening
- 2) Vaginal pressure measurement

# PLASMA (Centrifuge + Plasma Generator + Needle)





Bio Filler

# PRP Kit



# THIS BOX CONTAINS

- 3No 15ml Tubes
- 1No Pipette
- 1No Ampule
- 1No Blood Removal Cannula
- 1No Syringe (3 ml) with Needle



## Hair removal laser



#### SPECIFICATION:

Product name : Professional Diode Laser 755nm 808nm

1064nm Hair Removal Laser

Model : LHRPH

Wavelength : 755 + 808 + 1064nm

Laser Stack : USA Coherent 50 million shots guaranteed.

Energy Density : 1-300J/cm2
Pulse Width : 1-200ms
Repetition Rate : 1-10HZ

Skin Cooling Temperature: -5°C--5°C by cooled Sapphire

Spot Size : 15\*26mm 15\*35mm (Focus light maximumly)

Interface Language : English (Multi language)LCD Screen : 15.6" android system 4K


System Cooling Method : Water+radiator+fans+TEC condenser Input VAC : 220VAC/ 50-60Hz or 110VAC / 50-60Hz











- 1. 3 Wavelength: Treat all kind of hair color from black hair to white hair
- 2. Treat all skin types from white to dark skin
- 3. No pain and shorter treatment sessions
- 4. Effective and safe treatment for permanent hair removal

# Mesotherapy





Naked Device





Injecting Gun (Handle)

Device Components Kit







- Deap clean
- Lifting tightening
- re-hydration
- freckles fading skin rejuvenation
- noninvasive skin beauty

# SPECIFICATION:

Type : Mesotherapy Gun

Certification : E TUV ISO13485 FDA

Brand Name : DFBEAUTY / DFLASER, DFBEAUTY

Model Number : DF - MSG

Needle material : Titanium Alloy Treatment Head : 3size micro needle

Injection depth : 3.0mm

Technology : RF mesotheroy









# DERMA PEN WITH NEEDLE





Non-invasive Solution

# AR Magnetic Stimulation System

- Incontinence
- Sexual Dysfunction
- Chronic Pain
- Pelvic Organ Prolapse

#### **OVERVIEW**

- MEN & WOMEN & CHILDREN
- FULLY CLOTHED
- FREE OF OPERATION


#### **TECHNOLOGY**

FMS: Functional Magnetic Stimulation, the system generates pulsed magnetic field to cause pelvic floor muscles to performance thousands of effective contractions without direct contact with electrodes or probes.

#### **TECHNOLOGY**

SNMS: Sacral Nerve Magnetic Stimulation, the system sends mild electromagnetic stimulations on sacral roots to regulate pain and stress, improving urinary and bowel control.

# **ADVANTAGES**



- All Genders
   Benefits for Male, Female, and Children
- Adjustable Positions
   Allow effective stimulation
   in different area
- MEP Test
   Guide proper stimulation target

# Outstanding Performance

7 Tesla magnetic field

for deep penetration

100Hz pulse frequency

ensures effective treatment Liquid cooling system

for nonintermittent work

# All Genders

- Patients remain fully dressed, free of operation
- Predefined & customizable session plans for each patient
- Comfortable arm-chair design enables pleasant treatment experience



## KEGEL EXERCISER



A small comfortable probe is placed just a few centimetres into the vagina which sends gentle impulses to the correct pelvic floor muscles. This precise stimulation causes them to contract and get stronger quickly.

- Strengthens pelvic floor after childbirth
- Helps and stress incontinence naturally
- Improve sexual intimacy & satisfaction
- Discreet pelvic floor toning therphy
- Flattens tummy







Incontinence





mothers





sensitivity






Program



Exercise





INTIMATE AREA FILLER





# Glycopeel

Glycolic Acid

Available strengths 20%,35%,50%,70%

#### Indication:

- Aging skin Stimulate collagen production • Skin rejuvenation
- Fine lines Age spots
- Pigmentation Acne



# SALIC PEEL Salicylic Acid

Available strength 30%

#### Indication:

- Salicylic acid exfoliates skin
- Improve skin texture and color.
- Oily Skin

# MASTER PEEL

FOR PROFESSIONAL USE ONLY
Not for retail sale

# MASTER PEEL

Master Peel exfoliating treatment

TCA 30% & Phytic colour (blue).

#### Indication:

- Pigmentation Aging Skin
- Completely peels & removes dry flaky skin

100ml



# Ablaze Peel

Glycolic acid 50 % Lactic acid 10% Salicylic acid 4%

- Restores a radiant complexion
- Acts as epidermis stimulator
- Regulates oily skin
- Dilated pores
- Improves fine lines and wrinkles
- Even skin tone

50ml

#### **FUSION PEEL A**

Lactic 45% + Citric 5% + Pyruvic 5% + Arginine 5%



Wrinkles • Dark circle • Aging skin

#### **FUSION PEEL B**

Lactic 15% + Kojic 5% + Citric 5% + Sal 2%



 Skin lightening • Melasma Photodamaged skin

#### **FUSION PEEL C**

Sal 17% + Lactic 17% + Citric 8% + Gly 15%



 Anti Acne
 Anti Aging Freckles

**FUSION PEEL E** 

# **FUSION PEEL D**

Sal 20% + Man 10%



 Stretch Marks - Scarring Post - inflammatory acne pigmentation



50 ml

 Skin Rejuvenation - Superfacial Wrinkles









Lactic 20% Arginine 20%

#### Applications:

Reduces Wrinkles Increases Skin Hydration level Stimulate Collagen



Lactic 40%

#### Applications:

Treats Hyperpigmentation,
Dehydrated skin,
Premature aging,
Skin Lightening,
Lightens Sun Damaged
skin, Fades age spots,
Acne.



Lactic 80%

#### Applications:

Treating Dry, Rough &
Scaly skin, Excellent
Exfoliant and
Increases Skin
Hydration Level

50 ml

PELM1

50 ml

MANDLAZ - 20 Mandelic 20%



50 ml





Mandelic 20%

#### Applications:

Exfoliates Dead skin cells, Reduces the Apperance of Pores Prevent Breakout of Acne



Mandelic 30%

#### Applications:

50 ml

Simulates Collagen, Indicates for Oily and Seborrheic skin, Stimulates Collagen and Proteoglycan Synthesis.



Mandelic 40%

#### Applications:

Fine lines, Crystic acne, Large Pores, Darker Skin Tones and Sun Damaged skin, Loosen Blackheads.

50 ml

