

# 1470nm / 15 W Diode Laser with LVR and SUI Probe For Cosmetic Gynecology



SUI (Stress Urinary Incontinence) Handpiece

## THE BENEFITS OF LASER TREATMENT

#### A Simple Non Surgical Process that result in

- **Increased Collagen**
- Improved Muscle Tone
- **Improved Sensation**
- Reduced Vaginal Dryness
- **Enhanced Orgasms**
- OPD process without anesthesia
- No side effects
- Effective and painless

### **APPLICATIONS**

- Vaginal rejuvenation
- GSM-Genitourinary syndrome of menopause
- Vaginal atrophy
- Vaginal laxity
- Diseases associated to postpartum alteration
- HPV
- **Dryness**
- Itching
- Vaginismus

#### SPECIFICATION:

- Laser type
- Wavelength MAX Power
- Aiming Beam
- Operating Mode
- Pulse Time
- Beam delivery
- Optical Fiber Connector
- Emmition Initiation.
- Controller
- Display
- Cooling System

- : Diode, Semiconductor
- : 1470nm
- : 15 Watts
- : 635nm/Pulse, Single, Group
- : Continuous Pulse, Single, Group
- : 0.05ms-1000ms
- : Fiber Optic with SMA 905 Connector
- : SMA 905, accepts optical fibers having a core diameter from 200µm up to 600pm,
- : Footswitch
- : Microprocessor
- : 10.1" IPS with touch panel
- : Internal, air and thermoelectric cooling

- Power Supply of the Laser
- Power Supply of AC Adapter
- AC Adapter
- Laser Dimensions
- Laser weight
- Laser Case Dimensions
- Weight of Laser with Cases
- Class of Medical Device
- Laser Safety Class
- Electric Safety Class **Environmental Conditions During Work**

- : DC 24V/8.33A from the Separate AC
- Single phase 100-240 VAC: 50-60 Hz;
- DC 24v/8,33A Medical device class EN60601-1 Class |
- 27cm\* 11cm\*19 cm length\*width\*height
- 2.75kg
- $56\text{cm} \times 35\text{cm} \times 23\text{cm}$
- 09 kg
- IIB

- From +10 to 240C degree, relative humidity from 30% up to 60%





AMT, Plot No- A815, MIDC Industrial Area, Kopar Khairane, Navi Mumbai, Maharashtra 400709

**Office:** +91-9322673077 | +91-9892703077 **Office:** +91-8080321777 | +91-9321548914

Web: www.amtmedicals.com
Email: info@amtmedicals.com