



GYNECOLOGY

Equipment

COSMETIC GYNECOLOGY

TIZEN HIFU

High Intensity Focused Ultrasound



- 1) Sonic Vaginal rejuvenation
- 2) Sonic Vaginal tightening
- 2) Vaginal pressure measurement

How does Tizen HIFU works ?

The Tizen HIFU System uses a non-invasive ultrasonic focusing technique to target the mucosal lamia and muscle layer. Using ultrasonic waves as the energy source and taking advantage of its penetration and focality, a higher intensity of ultrasonic region, called focus region, is formed. In fraction of a second, the temperature of the region can reach slightly over 65 °C , so the collagen is reorganized and the normal tissue outside the focal region is undamaged. Therefore, the desired depth layer can obtain the ideal effect of collagen concentration, reorganization and regeneration.



Analyticals Medical
Technologies Pvt. Ltd.

AMT, Plot No- A815, MIDC Industrial Area,
Kopar Khairane, Navi Mumbai, Maharashtra 400709

Office : +91-9322673077 | +91-9892703077

Office : +91-8080321777 | +91-9321548914

Web : www.amtmedicals.com

Email : info@amtmedicals.com